



## Marco Compound # N1013

### 70 Durometer, Black, FDA Compliant Neoprene

### Technical Datasheet

#### **Common Names:**

Neoprene, Chloroprene, CR

#### **General Description:**

Neoprene was the first synthetic rubber developed commercially and exhibits generally good ozone, aging and chemical resistance. Neoprene rubbers contain Chlorine in the polymer to reduce the reactivity to many oxidizing agents, as well as to oil and flame. Neoprene have good abrasion and tear resistance and are suitable for use in heating and air conditioning systems (HVAC), refrigeration units and numerous dynamic applications. Please contact [engineering@marcorubber.com](mailto:engineering@marcorubber.com) for assistance in selecting a specialized compound when increased resistance to temperature, lubricants, or physical properties is required.

#### **Features:**

- FDA compliant
- Paraffin base mineral oil with low DPI, e.g. ASTM oil No. 1
- Silicone oil and grease
- Water and water solvents at low temperatures
- Refrigerants
- Ammonia
- Carbon dioxide
- Improved ozone, weathering and aging resistance compared with nitrile rubber.

#### **Limitations:**

- Aromatic hydrocarbons (benzene)
- Chlorinated hydrocarbons (trichloroethylene)
- Polar solvents (ketones, esters, ethers, acetones).

#### **Service Temperature:**

-40 to 250° F

#### **Specification:**

ASTM D2000 M3BC710 A14 B14 EO14 EO34 Z1 (Z1=FDA)

### PHYSICAL PROPERTY STANDARDS

ORIGINAL PROPERTIES	ASTM D2000 Requirements	Typical Test Results
Hardness, Shore A	70 +/- 5	72
Color	Black	Black
Tensile Strength, psi	1,450 min.	1,785
Ultimate Elongation, %	250 min.	320

Information within is believed to be accurate and reliable. However, Marco Rubber makes no warranty, expressed or implied, that parts supplied in this material will perform satisfactorily in specific applications. It's the customer's responsibility to evaluate prior to use.

<b>HEAT RESISTANCE – A14, ASTM D 573 (70 hrs. @ 100°C)</b>	<b>ASTM D2000 Requirements</b>	<b>Typical Test Results</b>
Hardness Change, points	+ 15	+8
Tensile Strength Change, %	- 15	+3
Ultimate Elongation Change, %	-40 max.	-21

<b>COMPRESSION SET – B14, ASTM D 325 Method B (22 hrs. @ 100°C)</b>	<b>ASTM D2000 Requirements</b>	<b>Typical Test Results</b>
Permanent Set %	35 max.	13

<b>FLUID RESISTANCE –ASTM #1 Oil – EO14, ASTM D 471 (70 hrs. @ 100°C)</b>	<b>ASTM D2000 Requirements</b>	<b>Typical Test Results</b>
Hardness Change, points	+/- 10	-9
Tensile Strength Change, %	-30 max.	-1
Ultimate Elongation Change, %	-30 max.	-21
Volume Change, %	-10 to +15	+1

<b>FLUID RESISTANCE – IRM 903 Oil, -EO34, ASTM D 471 (70 hrs. @ 100°C)</b>	<b>ASTM D2000 Requirements</b>	<b>Typical Test Results</b>
Tensile Strength Change, %	-60 max.	-57
Ultimate Elongation Change, %	-50 max.	-32
Volume Change, %	+100 max.	+58